

Product datasheet for **TA364590S**

Alpha 1 Acid Glycoprotein (ORM1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: Human testis and prostate tissue
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ORM1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	24 kDa
Gene Name:	orosomuroid 1
Database Link:	Entrez Gene 5004 Human P02763
Background:	This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.
Synonyms:	AGP-A; AGP1; ORM; Orosomuroid-1



[View online »](#)

Product images:



Gel: 12%SDS-PAGE

Lysate: 40 μ g

Lane 1-2: Human testis tissue

Human prostate tissue

Primary antibody: [TA364590] (ORM1 Antibody)
at dilution 1/133

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 2 minutes