

Product datasheet for **TA364690**

BMP15 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: A549 cells IHC: 50-200 Positive control: Human thyroid cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human BMP15
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	45 kDa
Gene Name:	bone morphogenetic protein 15
Database Link:	Entrez Gene 9210 Human O95972


Background: The protein encoded by this gene is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. The transforming growth factor-beta superfamily includes large families of growth and differentiation factors. It is thought that this protein may be involved in oocyte maturation and follicular development as a homodimer or by forming heterodimers with a related protein, Gdf9. Defects in this gene are the cause of ovarian dysgenesis 2.




[View online »](#)

Synonyms: BMP-15; GDF-9B; GDF9B; ODG2; POF4


Product images:




Gel: 10%SDS-PAGE
Lysate: 60 μ g
Lane: A549 cells
Primary antibody: TA364690 (BMP15 Antibody) at dilution 1/266.7
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 5 minutes




Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA364690 (BMP15 Antibody) at dilution 1/40 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA364690 (BMP15 Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA364690 (BMP15 Antibody) at dilution 1/40 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA364690 (BMP15 Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: $\times 200$)