

Product datasheet for **TA364756S**


DHRS3 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Mouse liver tissue and A375 cells IHC: 50-200 Positive control: Human breast cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human DHRS3
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	34 kDa
Gene Name:	dehydrogenase/reductase 3
Database Link:	Entrez Gene 9249 Human O75911
Background:	Short-chain dehydrogenases/reductases (SDRs), such as DHRS3, catalyze the oxidation/reduction of a wide range of substrates, including retinoids and steroids.
Synonyms:	DD83.1; RDH17; retSDR1; Rsd1; SDR1; SDR16C1




[View online »](#)

Product images:


Gel: 10%SDS-PAGE
Lysate: 40 µg
Lane 1-2: Mouse liver tissue
A375 cells
Primary antibody: [TA364756] (DHRS3 Antibody)
at dilution 1/550
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 20 seconds




Immunohistochemistry of paraffin-embedded
Human breast cancer tissue using [TA364756]
(DHRS3 Antibody) at dilution 1/30 (Original
magnification: ×200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA364756] (DHRS3 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using [TA364756] (DHRS3 Antibody) at dilution 1/30 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using [TA364756] (DHRS3 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: $\times 200$)