


## Product datasheet for **TA364885S**

### MFAP3L Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-300 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human MFAP3L
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	microfibrillar associated protein 3 like
Database Link:	<a href="#">Entrez Gene 9848 Human O75121</a>
Background:	MFAP3L (microfibrillar-associated protein 3-like), also known as HSD39 or testis development protein NYD-SP9, is a 409 amino acid single-pass type I cell membrane protein that contains one Ig-like (immunoglobulin-like) domain. Found primarily in testis, MFAP3L is encoded by a gene that is located on chromosome 4 and is expressed as three isoforms due to alternative splicing events. Representing approximately 6% of the human genome, chromosome 4 contains nearly 900 genes, one of which is the Huntingtin gene, which is found to encode an expanded glutamine tract in cases of Huntington's disease.
Synonyms:	KIAA0626; NYD-sp9


[View online »](#)

**Product images:**



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA364885] (MFAP3L Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA364885] (MFAP3L Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA364885] (MFAP3L Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA364885] (MFAP3L Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)