

Product datasheet for **TA364895S**


NADE (BEX3) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 200-400 Positive control: Human breast cancer Predicted cell location: Cytoplasm or Nucleus
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human BEX3
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	brain expressed X-linked 3
Database Link:	Entrez Gene 27018 Human Q00994
Background:	May be a signaling adapter molecule involved in p75 ^{NTR} -mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death (By similarity). May play an important role in the pathogenesis of neurogenetic diseases.




[View online »](#)

Product images:


Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA364895] (BEX3 Antibody) at dilution 1/190 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA364895] (BEX3 Antibody) at dilution 1/190, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA364895] (BEX3 Antibody) at dilution 1/190 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA364895] (BEX3 Antibody) at dilution 1/190, treated with fusion protein. (Original magnification: $\times 200$)