

Product datasheet for **TA365240**


RUFY4 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human lung cancer Predicted cell location: Cytoplasm or Nucleus
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human RUFY4
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	RUN and FYVE domain containing 4
Database Link:	Entrez Gene 285180 Human Q6ZNE9
Background:	The function of this protein remains unknown.
Synonyms:	FLJ46536




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA365240 (RUFY4 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA365240 (RUFY4 Antibody) at dilution 1/35, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA365240 (RUFY4 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA365240 (RUFY4 Antibody) at dilution 1/35, treated with fusion protein. (Original magnification: ×200)