

Product datasheet for **TA365245**

Marapsin2 (PRSS38) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Secreted
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	protease, serine 38
Database Link:	Entrez Gene 339501 Human A1L453
Background:	The function of this protein remains unknown.
Synonyms:	MPN2




[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA365245 (PRSS38 Antibody) at dilution 1/30 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA365245 (PRSS38 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA365245 (PRSS38 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA365245 (PRSS38 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)