

## Product datasheet for **TA365300**


### DYDC2 Rabbit Polyclonal Antibody

#### Product data:


|                       |                                                                                                                                                                                                                             |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                                                                                                                                                          |
| Applications:         | IHC                                                                                                                                                                                                                         |
| Recommended Dilution: | IHC: 25-100<br>Positive control: Human brain<br>Predicted cell location: Cytoplasm and Cell membrane                                                                                                                        |
| Reactivity:           | Human                                                                                                                                                                                                                       |
| Host:                 | Rabbit                                                                                                                                                                                                                      |
| Isotype:              | IgG                                                                                                                                                                                                                         |
| Clonality:            | Polyclonal                                                                                                                                                                                                                  |
| Immunogen:            | Full length fusion protein                                                                                                                                                                                                  |
| Formulation:          | pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol                                                                                                                                                                            |
| Concentration:        | lot specific                                                                                                                                                                                                                |
| Purification:         | Antigen affinity purification                                                                                                                                                                                               |
| Conjugation:          | Unconjugated                                                                                                                                                                                                                |
| Storage:              | Store at -20°C.                                                                                                                                                                                                             |
| Stability:            | 1 year                                                                                                                                                                                                                      |
| Gene Name:            | DPY30 domain containing 2                                                                                                                                                                                                   |
| Database Link:        | <a href="#">Entrez Gene 84332 Human Q96IM9</a>                                                                                                                                                                              |
| Background:           | This gene encodes a member of a family of proteins that contains a DPY30 domain. This gene locus overlaps with a closely related gene on the opposite strand. Alternative splicing results in multiple transcript variants. |
| Synonyms:             | bA36D19.6; MGC16186; OTTHUMP00000019966; OTTHUMP00000019967                                                                                                                                                                 |



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human brain tissue using TA365300 (DYDC2 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using TA365300 (DYDC2 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)