

Product datasheet for **TA365312S**

TESC Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human prostate cancer Predicted cell location: Cytoplasm and Cell membrane
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human TESC
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	tescalcin
Database Link:	Entrez Gene 54997 Human Q96BS2




[View online »](#)

Background:


Tescalcin, also known as TESC, TSC or CHP3, is a 267 amino acid protein that contains one EF-hand domain and is expressed abundantly in adult heart tissue. Using calcium as a cofactor, Tescalcin interacts with NHE-1 and functions to couple the activation of the ERK cascade with the expression of Ets proteins during megakaryocytic differentiation. Human Tescalcin shares 97% sequence identity with its mouse counterpart, suggesting a conserved role between species. Multiple isoforms of Tescalcin exist due to alternative splicing events. Functions as an integral cofactor in cell pH regulation by controlling plasma membrane-type Na⁺/H⁺ exchange activity. Promotes the maturation, transport, cell surface stability and exchange activity of SLC9A1/NHE1 at the plasma membrane. Promotes the induction of hematopoietic stem cell differentiation toward megakaryocytic lineage. Essential for the coupling of ERK cascade activation with the expression of ETS family genes in megakaryocytic differentiation. Also involved in granulocytic differentiation in a ERK-dependent manner. Inhibits the phosphatase activity of calcineurin.

Synonyms:


CHP3; FLJ20607; OTTHUMP00000185561; tescalcin; TSC

Product images:


Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA365312] (TESC Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA365312] (TESC Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA365312] (TESC Antibody) at dilution 1/25 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA365312] (TESC Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: $\times 200$)