

Product datasheet for **TA365331**


ZNF22 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: Jurkat cell and Human fetal brain tissue lysates IHC: 25-100 Positive control: Human prostate cancer Predicted cell location: Cytoplasm and Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZNF22
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	26 kDa
Gene Name:	zinc finger protein 22
Database Link:	Entrez Gene 7570 Human P17026
Background:	Binds DNA through the consensus sequence 5'-CAATG-3'. May be involved in transcriptional regulation and may play a role in tooth formation (By similarity).
Synonyms:	HKR-T1; KOX15; KROX26; Zfp422; ZNF422



[View online »](#)

Product images:

Gel: 12%SDS-PAGE


Lysate: 40 μ g

Lane 1-2: Jurkat cell and Human fetal brain tissue lysates


Primary antibody: TA365331 (ZNF22 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution


Exposure time: 20 seconds




Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using TA365331 (ZNF22 Antibody) at dilution 1/25 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using TA365331 (ZNF22 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA365331 (ZNF22 Antibody) at dilution 1/25 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA365331 (ZNF22 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: $\times 200$)