

Product datasheet for **TA365378S**


ZKSCAN3 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-300 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm, Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZKSCAN3
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	zinc finger with KRAB and SCAN domains 3
Database Link:	Entrez Gene 80317 Human Q9BRR0
Background:	Transcriptional factor that binds to the consensus sequence 5'-[GT][AG][AGT]GGGG-3' and acts as a repressor of autophagy. Specifically represses expression of genes involved in autophagy and lysosome biogenesis/function such as MAP1LC3B, ULK1 or WIPI2. Associates with chromatin at the ITGB4 and VEGF promoters. Also acts as a transcription activator and promotes cancer cell progression and/or migration in various tumors and myelomas.
Synonyms:	dj874C20.1; FLJ33906; KIAA0426; OTTHUMP00000016204; ZF47; Zfp-47; Zfp47; ZFP306; ZNF306; ZNF309; ZSCAN13



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA365378] (ZKSCAN3 Antibody) at dilution 1/125 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA365378] (ZKSCAN3 Antibody) at dilution 1/125, treated with fusion protein. (Original magnification: $\times 200$)