

Product datasheet for **TA365379S**

ZMAT2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: Jurkat,231,Hela,k562 and A172 cell,Mouse brain tissue lysates IHC: 30-150 Positive control: Human tonsil Predicted cell location: Cytoplasm and Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	24 kDa
Gene Name:	zinc finger matrin-type 2
Database Link:	Entrez Gene 153527 Human Q96NC0
Background:	The function of ZMAT2 remains unknown.
Synonyms:	FLJ31121




[View online »](#)


Product images:




Gel: 12%SDS-PAGE
Lysate: 40 μ g
Lane 1-6: Jurkat
231
Hela
k562 and A172 cell
Mouse brain tissue lysates
Primary antibody: [TA365379] (ZMAT2 Antibody)
at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 3 seconds




Immunohistochemistry of paraffin-embedded
Human tonsil tissue using [TA365379] (ZMAT2
Antibody) at dilution 1/40 (Original magnification:
 \times 200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA365379] (ZMAT2 Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA365379] (ZMAT2 Antibody) at dilution 1/40 (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA365379] (ZMAT2 Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: x200)