

Product datasheet for **TA365406S**


ZNF690 (ZSCAN29) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human prostate cancer Predicted cell location: Nucleus
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZSCAN29
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	zinc finger and SCAN domain containing 29
Database Link:	Entrez Gene 146050 Human Q8IWY8
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ35867; MGC129894; MGC129895; Zfp690; ZNF690

[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA365406] (ZSCAN29 Antibody) at dilution 1/30 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA365406] (ZSCAN29 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA365406] (ZSCAN29 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA365406] (ZSCAN29 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)