

Product datasheet for **TA365516**

BCO2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human brain Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human BCO2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	beta-carotene oxygenase 2
Database Link:	Entrez Gene 83875 Human Q9BYV7
Background:	This gene encodes an enzyme which oxidizes carotenoids such as beta-carotene during the biosynthesis of vitamin A. Multiple transcript variants encoding different isoforms have been found for this gene.
Synonyms:	B-DIOX-II; BCDO2; FLJ34464; OTTHUMP00000196941




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human brain tissue using TA365516 (BCO2 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using TA365516 (BCO2 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)