

## Product datasheet for **TA365785S**

### UBXN1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: 293T□HEPG2□Hela□LOVO and 231 cell□Mouse brain tissue lysates IHC: 25-100 Positive control: Human thyroid cancer Predicted cell location: Cytoplasm and Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	33 kDa
Gene Name:	UBX domain protein 1
Database Link:	<a href="#">Entrez Gene 51035 Human Q04323</a>
Background:	Ubiquitin-binding protein that interacts with the BRCA1-BARD1 heterodimer, and regulates its activity. Specifically binds 'Lys-6'-linked polyubiquitin chains. Interaction with autoubiquitinated BRCA1, leads to inhibit the E3 ubiquitin-protein ligase activity of the BRCA1-BARD1 heterodimer. Component of a complex required to couple deglycosylation and proteasome-mediated degradation of misfolded proteins in the endoplasmic reticulum that are retrotranslocated in the cytosol.
Synonyms:	2B28; LOC51035; SAKS1; UBXD10




[View online »](#)


## Product images:




Gel: 12%SDS-PAGE  
Lysate: 40  $\mu$ g  
Lane 1-6: 293T  
HEPG2  
Hela  
LOVO and 231 cell  
Mouse brain tissue lysates  
Primary antibody: [TA365785] (UBXN1 Antibody)  
at dilution 1/250  
Secondary antibody: Goat anti rabbit IgG at  
1/8000 dilution  
Exposure time: 10 seconds




Immunohistochemistry of paraffin-embedded  
Human thyroid cancer tissue using [TA365785]  
(UBXN1 Antibody) at dilution 1/20 (Original  
magnification:  $\times$ 200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA365785] (UBXN1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA365785] (UBXN1 Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA365785] (UBXN1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: ×200)