

Product datasheet for **TA365827S**


RMND5B Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human esophagus cancer Predicted cell location: Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human RMND5B
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	required for meiotic nuclear division 5 homolog B
Database Link:	Entrez Gene 64777 Human Q96G75
Background:	The function of this protein remains unknown.
Synonyms:	DKFZp434K0926; FLJ22318



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA365827] (RMND5B Antibody) at dilution 1/50 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA365827] (RMND5B Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: $\times 200$)