

## Product datasheet for **TA365870S**


### PCTAIRE3 (CDK18) Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 30-150 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human CDK18
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	cyclin-dependent kinase 18
Database Link:	<a href="#">Entrez Gene 5129 Human Q07002</a>
Background:	May play a role in signal transduction cascades in terminally differentiated cells.
Synonyms:	PCTAIRE; PCTAIRE3; PCTK3



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA365870] (CDK18 Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA365870] (CDK18 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: ×200)