

Product datasheet for **TA365939**


PDRG1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human PDRG1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	p53 and DNA damage regulated 1
Database Link:	Entrez Gene 81572 Human Q9NUG6
Background:	May play a role in chaperone-mediated protein folding.
Synonyms:	C20orf126; dj310013.3; PDRG



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA365939 (PDRG1 Antibody) at dilution 1/80 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA365939 (PDRG1 Antibody) at dilution 1/80, treated with fusion protein. (Original magnification: ×200)