

Product datasheet for **TA365986S**

NIPSNAP3A Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human cervical cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	nipsnap homolog 3A
Database Link:	Entrez Gene 25934 Human Q9UFN0
Background:	NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport (Buechler et al., 2004 [PubMed 15177564]).
Synonyms:	DKFZp564D177; FLJ13953; HSPC299; MGC14553; NIPSNAP4; TASSC




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA365986] (NIPSNAP3A Antibody) at dilution 1/40 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA365986] (NIPSNAP3A Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: ×200)