

Product datasheet for **TA366115S**

NPC2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Mouse kidney tissue lysate IHC: 50-200 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human NPC2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	17 kDa
Gene Name:	NPC intracellular cholesterol transporter 2
Database Link:	Entrez Gene 10577 Human P61916
Background:	This gene encodes a protein containing a lipid recognition domain. The encoded protein may function in regulating the transport of cholesterol through the late endosomal/lysosomal system. Mutations in this gene have been associated with Niemann-Pick disease, type C2 and frontal lobe atrophy.
Synonyms:	EDDM1; HE1; MGC1333; NP-C2




[View online »](#)


Product images:




Gel: 12%SDS-PAGE
Lysate: 40 μ g
Lane: Mouse kidney tissue lysate
Primary antibody: [TA366115] (NPC2 Antibody) at dilution 1/750
Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution
Exposure time: 3 minutes




Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA366115] (NPC2 Antibody) at dilution 1/50 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA366115] (NPC2 Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA366115] (NPC2 Antibody) at dilution 1/50 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA366115] (NPC2 Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: $\times 200$)