

Product datasheet for **TA366189S**

ZNF688 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human muscle tissue lysate IHC: 50-200 Positive control: Human prostate cancer Predicted cell location: Nucleus
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZNF688
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	31 kDa
Gene Name:	zinc finger protein 688
Database Link:	Entrez Gene 146542 Human P0C7X2
Background:	May be involved in transcriptional regulation.
Synonyms:	ZNF688




[View online »](#)


Product images:



Gel: 12%SDS-PAGE
Lysate: 40 μ g
Lane: Human muscle tissue lysate
Primary antibody: [TA366189] (ZNF688 Antibody) at dilution 1/800
Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution
Exposure time: 3 minutes



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA366189] (ZNF688 Antibody) at dilution 1/110 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA366189] (ZNF688 Antibody) at dilution 1/110, treated with fusion protein. (Original magnification: ×200)