

Product datasheet for **TA366254S**


TENT5C Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-300 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human TENT5C
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	family with sequence similarity 46 member C
Database Link:	Entrez Gene 54855 Human Q5VWP2
Background:	Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. (Microbial infection) Seems to enhance replication of some viruses, including yellow fever virus, in response to type I interferon.
Synonyms:	FLJ20202



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366254] (TENT5C Antibody) at dilution 1/110 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366254] (TENT5C Antibody) at dilution 1/110, treated with fusion protein. (Original magnification: ×200)