

Product datasheet for **TA366277S**

KIR3DL1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1000-5000 WB positive control: Human fetal liver tissue and Human liver tissue lysates IHC: 50-300 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human KIR3DL1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	49 kDa
Gene Name:	killer cell immunoglobulin like receptor, three Ig domains and long cytoplasmic tail 1
Database Link:	Entrez Gene 3811 Human P43629




[View online »](#)

Background:

Killer cell immunoglobulin-like receptors (KIRs) are transmembrane glycoproteins expressed by natural killer cells and subsets of T cells. The KIR genes are polymorphic and highly homologous and they are found in a cluster on chromosome 19q13.4 within the 1 Mb leukocyte receptor complex (LRC). The gene content of the KIR gene cluster varies among haplotypes, although several framework genes are found in all haplotypes (KIR3DL3, KIR3DP1, KIR3DL4, KIR3DL2). The KIR proteins are classified by the number of extracellular immunoglobulin domains (2D or 3D) and by whether they have a long (L) or short (S) cytoplasmic domain. KIR proteins with the long cytoplasmic domain transduce inhibitory signals upon ligand binding via an immune tyrosine-based inhibitory motif (ITIM) while KIR proteins with the short cytoplasmic domain lack the ITIM motif and instead associate with the TYRO protein tyrosine kinase binding protein to transduce activating signals. The ligands for several KIR proteins are subsets of HLA class I molecules; thus KIR proteins are thought to play an important role in regulation of the immune response.

Synonyms:

AMB11; CD158E; CD158E1; CD158e1/2; CD158e2; cl-2; cl-11; KIR; KIR3DS1; MGC119726; MGC119728; MGC126589; MGC126591; NK-receptor; NKAT-3; NKAT3; NKB1; NKB1B

Product images:

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-2: Human fetal liver tissue and Human liver tissue lysates


Primary antibody: [TA366277] (KIR3DL1 Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution


Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366277] (KIR3DL1 Antibody) at dilution 1/105 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366277] (KIR3DL1 Antibody) at dilution 1/105, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA366277] (KIR3DL1 Antibody) at dilution 1/105 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA366277] (KIR3DL1 Antibody) at dilution 1/105, treated with fusion protein. (Original magnification: $\times 200$)