


Product datasheet for **TA366294S**

ZNF256 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50–300 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZNF256
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at –20°C.
Stability:	1 year
Gene Name:	zinc finger protein 256
Database Link:	Entrez Gene 10172 Human Q9Y2P7
Background:	Transcriptional repressor that plays a role in cell proliferation. Requires TRIM28 for its activity.
Synonyms:	BMZF-3; BMZF3




Product images:


Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA366294] (ZNF256 Antibody) at dilution 1/70 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA366294] (ZNF256 Antibody) at dilution 1/70, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366294] (ZNF256 Antibody) at dilution 1/70 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366294] (ZNF256 Antibody) at dilution 1/70, treated with fusion protein. (Original magnification: ×200)