

Product datasheet for **TA366305**

ZNFN1A2 (IKZF2) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1000-5000 WB positive control: Human skin tissue, A172 cell, Human bladder transitional cell carcinoma grade 2-3 tissue lysates
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human IKZF2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	60 kDa
Gene Name:	IKAROS family zinc finger 2
Database Link:	Entrez Gene 22807 Human Q9UKS7
Background:	This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This protein forms homo- or heterodimers with other Ikaros family members, and is thought to function predominantly in early hematopoietic development. Multiple transcript variants encoding different isoforms have been found for this gene, but the biological validity of some variants has not been determined.
Synonyms:	helios; MGC34330; OTTHUMP00000163964; ZNF1A2; ZNFN1A2



[View online »](#)

Product images:



Gel: 8%SDS-PAGE

Lysate: 40 μ g

Lane: Human skin tissue

A172 cell

Human bladder transitional cell carcinoma grade 2-3 tissue lysates

Primary antibody: TA366305 (IKZF2 Antibody) at dilution 1/950

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 1 minute