

Product datasheet for **TA366365S**

MRPL42 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: 293T, Hela, HepG2 and Jurkat cell lysates IHC: 50-200 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human MRPL42
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	17 kDa
Gene Name:	mitochondrial ribosomal protein L42
Database Link:	Entrez Gene 28977 Human Q9Y6G3




[View online »](#)

Background:


Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a protein identified as belonging to both the 28S and the 39S subunits. Alternative splicing results in multiple transcript variants. Pseudogenes corresponding to this gene are found on chromosomes 4q, 6p, 6q, 7p, and 15q.

Synonyms:


HSPC204; L42mt; MRP-L31; MRP-L42; MRP-S32; MRPL31; MRPS32; PTD007; RPML31; S32mt

Product images:


Gel: 12%SDS-PAGE
Lysate: 40 µg
Lane 1-4: 293T
HeLa
HepG2 and Jurkat cell lysates
Primary antibody: [TA366365] (MRPL42 Antibody)
at dilution 1/650
Secondary antibody: Goat anti rabbit IgG at
1/5000 dilution
Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded
Human liver cancer tissue using [TA366365]
(MRPL42 Antibody) at dilution 1/60 (Original
magnification: x200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366365] (MRPL42 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA366365] (MRPL42 Antibody) at dilution 1/60 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA366365] (MRPL42 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: $\times 200$)