

## Product datasheet for **TA366468S**


### **C11orf54 Rabbit Polyclonal Antibody**

#### **Product data:**


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-200 Positive control: Human liver cancer Predicted cell location: Cytoplasm and Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human C11orf54
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	chromosome 11 open reading frame 54
Database Link:	<a href="#">Entrez Gene 28970 Human Q9H0W9</a>
Background:	Exhibits ester hydrolase activity on the substrate p-nitrophenyl acetate.
Synonyms:	LP4947; PTD012



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366468] (C11orf54 Antibody) at dilution 1/70 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA366468] (C11orf54 Antibody) at dilution 1/70, treated with fusion protein. (Original magnification:  $\times 200$ )