

## Product datasheet for **TA366494**


### **DHRS7B Rabbit Polyclonal Antibody**

#### **Product data:**


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human DHRS7B
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	dehydrogenase/reductase 7B
Database Link:	<a href="#">Entrez Gene 25979 Human Q6IAN0</a>
Background:	This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.
Synonyms:	CGI-93; DKFZp566O084; MGC8916; SDR32C1



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA366494 (DHR57B Antibody) at dilution 1/80 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA366494 (DHR57B Antibody) at dilution 1/80, treated with fusion protein. (Original magnification:  $\times 200$ )