

Product datasheet for **TA366509S**

HEATR1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-100 Positive control: Human prostate cancer Predicted cell location: Nucleus
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human HEATR1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	HEAT repeat containing 1
Database Link:	Entrez Gene 55127 Human Q9H583
Background:	Ribosome biogenesis factor. Involved in nucleolar processing of pre-18S ribosomal RNA. Required for optimal pre-ribosomal RNA transcription by RNA polymerase I.
Synonyms:	BAP28; FLJ10359; MGC72083




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA366509] (HEATR1 Antibody) at dilution 1/80 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA366509] (HEATR1 Antibody) at dilution 1/80, treated with fusion protein. (Original magnification: $\times 200$)