

Product datasheet for **TA366551S**


DEPTOR Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1000-5000 WB positive control: SKOV3 cell lysate IHC: 150-300 Positive control: Human prostate cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human DEPTOR
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	46 kDa
Gene Name:	DEP domain containing MTOR-interacting protein
Database Link:	Entrez Gene 64798 Human Q8TB45
Background:	Negative regulator of the mTORC1 and mTORC2 signaling pathways. Inhibits the kinase activity of both complexes.
Synonyms:	DEP.6; DEPDC6




[View online »](#)

Product images:

Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane: SKOV3 cell lysate
Primary antibody: [TA366551] (DEPTOR Antibody)
at dilution 1/1000
Secondary antibody: Goat anti rabbit IgG at
1/5000 dilution
Exposure time: 1 minute



Immunohistochemistry of paraffin-embedded
Human prostate cancer tissue using [TA366551]
(DEPTOR Antibody) at dilution 1/160 (Original
magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA366551] (DEPTOR Antibody) at dilution 1/160, treated with fusion protein. (Original magnification: $\times 200$)