

Product datasheet for **TA367296S**


STK32C Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human STK32C
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	serine/threonine kinase 32C
Database Link:	Entrez Gene 282974 Human Q86UX6
Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.
Synonyms:	MGC23665; PKE; RP11-140A10.1; YANK3



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA367296] (STK32C Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA367296] (STK32C Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)