

## Product datasheet for **TA367568S**

### NUDT8 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human NUDT8
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	nudix hydrolase 8
Database Link:	<a href="#">Entrez Gene 254552 Human Q8WV74</a>
Background:	Probably mediates the hydrolysis of some nucleoside diphosphate derivatives.
Synonyms:	FLJ41567




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA367568] (NUDT8 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA367568] (NUDT8 Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)