

Product datasheet for **TA367709**

LRRC49 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-50 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human LRRC49
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	leucine rich repeat containing 49
Database Link:	Entrez Gene 54839 Human Q8IUZ0
Background:	The function of this protein remains unknown.
Synonyms:	FLJ20156; PGs4




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA367709 (LRRC49 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA367709 (LRRC49 Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: ×200)