

Product datasheet for **TA367905**

C16orf88 (KNOP1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human esophagus cancer Predicted cell location: Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human KNOP1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	lysine-rich nucleolar protein 1
Database Link:	Entrez Gene 400506 Human Q1ED39
Background:	poly(A) RNA binding
Synonyms:	101F10.1; C16orf88; FAM191A; TSG118




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA367905 (KNOP1 Antibody) at dilution 1/85 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA367905 (KNOP1 Antibody) at dilution 1/85, treated with synthetic peptide. (Original magnification: $\times 200$)