

Product datasheet for **TA367930S**


Calpain 14 (CAPN14) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Jurkat and A549 cell lysates IHC: 50-100 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human CAPN14
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	79 kDa
Gene Name:	calpain 14
Database Link:	Entrez Gene 440854 Human A8MX76
Background:	Calpains are a family of cytosolic calcium-activated cysteine proteases involved in a variety of cellular processes including apoptosis, cell division, modulation of integrin-cytoskeletal interactions, and synaptic plasticity. CAPN14 belongs to the calpain large subunit family.
Synonyms:	FLJ54788

[View online »](#)


Product images:




Gel: 8%SDS-PAGE
Lysate: 40 µg
Lane 1-2: Jurkat and A549 cell lysates
Primary antibody: [TA367930] (CAPN14 Antibody) at dilution 1/1000
Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution
Exposure time: 5 minutes




Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA367930] (CAPN14 Antibody) at dilution 1/90 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA367930] (CAPN14 Antibody) at dilution 1/90, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA367930] (CAPN14 Antibody) at dilution 1/90 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA367930] (CAPN14 Antibody) at dilution 1/90, treated with synthetic peptide. (Original magnification: ×200)