

Product datasheet for **TA368070S**


NCBP1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: Jurkat□NIH/3T3□293T□Hela□HT-29 and HEPG2 cell lysates IHC: 25-100 Positive control: Human lung cancer Predicted cell location: Nucleus or Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human NCBP1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	92 kDa
Gene Name:	nuclear cap binding protein subunit 1
Database Link:	Entrez Gene 4686 Human Q09161
Background:	The product of this gene is a component of the nuclear cap-binding protein complex (CBC), which binds to the monomethylated 5' cap of nascent pre-mRNA in the nucleoplasm. The encoded protein promotes high-affinity mRNA-cap binding and associates with the CTD of RNA polymerase II. The CBC promotes pre-mRNA splicing, 3'-end processing, RNA nuclear export, and nonsense-mediated mRNA decay.
Synonyms:	CBP80; MGC2087; NCBP



[View online »](#)

Product images:

Gel: 6%SDS-PAGE

Lysate: 40 µg

Lane 1-6: Jurkat

NIH/3T3

293T


HeLa

HT-29 and HEPG2 cell lysates


Primary antibody: [TA368070] (NCBP1 Antibody)
at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution


Exposure time: 40 seconds




Immunohistochemistry of paraffin-embedded
Human lung cancer tissue using [TA368070]
(NCBP1 Antibody) at dilution 1/20 (Original
magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using [TA368070] (NCBP1 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA368070] (NCBP1 Antibody) at dilution 1/20 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA368070] (NCBP1 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: $\times 200$)