

Product datasheet for **TA368146**

CYLN2 (CLIP2) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: NIH/3T3 and LOVO cell lysates IHC: 20-100 Positive control: Human thyroid cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human CLIP2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	116 kDa
Gene Name:	CAP-Gly domain containing linker protein 2
Database Link:	Entrez Gene 7461 Human Q9UDT6

Background: The protein encoded by this gene belongs to the family of cytoplasmic linker proteins, which have been proposed to mediate the interaction between specific membranous organelles and microtubules. This protein was found to associate with both microtubules and an organelle called the dendritic lamellar body. This gene is hemizygotously deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. Alternative splicing of this gene generates 2 transcript variants.


[View online »](#)

Synonyms:

CLIP; CLIP-115; CYLN2; KIAA0291; MGC11333; OTTHUMP00000160724; WBSCR3; WBSCR4; WSCR3; WSCR4

Product images:


Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-2: NIH/3T3 and LOVO cell lysates


Primary antibody: TA368146 (CLIP2 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution


Exposure time: 40 seconds




Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA368146 (CLIP2 Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA368146 (CLIP2 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA368146 (CLIP2 Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA368146 (CLIP2 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: ×200)