

Product datasheet for **TA368195**


DNAJC9 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: Human tonsil tissue lysate IHC: 20-100 Positive control: Human colorectal cancer Predicted cell location: Nucleus or Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human DNAJC9
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	30 kDa
Gene Name:	DnaJ heat shock protein family (Hsp40) member C9
Database Link:	Entrez Gene 23234 Human Q8WXX5
Background:	May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8.
Synonyms:	HDJC9; JDD1; KIAA0974; SB73




[View online »](#)

Product images:


Gel: 12%SDS-PAGE
Lysate: 40 μ g
Lane: Human tonsil tissue lysate
Primary antibody: TA368195 (DNAJC9 Antibody)
at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded
Human colorectal cancer tissue using TA368195
(DNAJC9 Antibody) at dilution 1/20 (Original
magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA368195 (DNAJC9 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA368195 (DNAJC9 Antibody) at dilution 1/20 (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA368195 (DNAJC9 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: x200)