

## Product datasheet for **TA368279S**


### **EHMT2/G9A (EHMT2) Rabbit Polyclonal Antibody**

#### **Product data:**


<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	IHC
<b>Recommended Dilution:</b>	IHC: 25-100 Positive control: Human ovarian cancer Predicted cell location: Nucleus
<b>Reactivity:</b>	Human, Mouse
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	Synthetic peptide of human EHMT2
<b>Formulation:</b>	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification:</b>	Antigen affinity purification
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	1 year
<b>Gene Name:</b>	euchromatic histone lysine methyltransferase 2
<b>Database Link:</b>	<a href="#">Entrez Gene 10919 Human Q96KQ7</a>
<b>Background:</b>	This gene encodes a methyltransferase that methylates lysine residues of histone H3. Methylation of H3 at lysine 9 by this protein results in recruitment of additional epigenetic regulators and repression of transcription. This gene was initially thought to be two different genes, NG36 and G9a, adjacent to each other in the HLA locus. Alternative splicing results in multiple transcript variants.
<b>Synonyms:</b>	BAT8; C6orf30; DKFZp686H08213; Em:AF134726.3; FLJ35547; G9A; KMT1C; NG36; NG36/G9a




[View online »](#)

**Product images:**


Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using [TA368279] (EHMT2 Antibody) at dilution 1/20 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using [TA368279] (EHMT2 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA368279] (EHMT2 Antibody) at dilution 1/20 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA368279] (EHMT2 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification:  $\times 200$ )