

## Product datasheet for **TA368490**

### Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNPA1L2) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 500-2000 WB positive control: Jurkat and 293T cell lysates
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human HNRNPA1L2
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	34 kDa
Gene Name:	heterogeneous nuclear ribonucleoprotein A1-like 2
Database Link:	<a href="#">Entrez Gene 144983 Human</a> <a href="#">Q32P51</a>
Background:	Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.
Synonyms:	HNRNPA1L; LOC144983; MGC102957; OTTHUMP00000018460



[View online »](#)

## Product images:



Gel: 12%SDS-PAGE

Lysate: 40  $\mu$ g

Lane 1-2: Jurkat and 293T cell lysates

Primary antibody: TA368490 (HNRNPA1L2 Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 1 minute