

Product datasheet for **TA368499S**

HLA DOB (HLA-DOB) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: HepG2, Hela, 293T cell lysate IHC: 30-150 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human HLA-DOB
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	31 kDa
Gene Name:	major histocompatibility complex, class II, DO beta
Database Link:	Entrez Gene 3112 Human P13765




[View online »](#)

Background:


HLA-DOB belongs to the HLA class II beta chain paralogues. This class II molecule is a heterodimer consisting of an alpha (DOA) and a beta chain (DOB), both anchored in the membrane. It is located in intracellular vesicles. DO suppresses peptide loading of MHC class II molecules by inhibiting HLA-DM. Class II molecules are expressed in antigen presenting cells (APC: B lymphocytes, dendritic cells, macrophages). The beta chain is approximately 26-28 kDa and its gene contains 6 exons. Exon one encodes the leader peptide, exons 2 and 3 encode the two extracellular domains, exon 4 encodes the transmembrane domain and exon 5 encodes the cytoplasmic tail.

Product images:


Gel: 12%SDS-PAGE
Lysate: 40 µg
Lane 1-3: HepG2
Hela
293T cell lysate
Primary antibody: [TA368499] (HLA-DOB
Antibody) at dilution 1/400
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded
Human tonsil tissue using [TA368499] (HLA-DOB
Antibody) at dilution 1/40 (Original magnification:
×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA368499] (HLA-DOB Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA368499] (HLA-DOB Antibody) at dilution 1/40 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA368499] (HLA-DOB Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: $\times 200$)