

## Product datasheet for **TA368552S**


### INSM2 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Mouse heart tissue lysate IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human INSM2
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	59 kDa
Gene Name:	INSM transcriptional repressor 2
Database Link:	<a href="#">Entrez Gene 84684 Human Q96T92</a>
Background:	May function as a growth suppressor or tumor suppressor in liver cells and in certain neurons.
Synonyms:	IA-6; IA6; mlt1




[View online »](#)

**Product images:**


Gel: 8%SDS-PAGE  
Lysate: 40 µg  
Lane: Mouse heart tissue lysate  
Primary antibody: [TA368552] (INSM2 Antibody) at dilution 1/400  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA368552] (INSM2 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA368552] (INSM2 Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA368552] (INSM2 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA368552] (INSM2 Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: ×200)