

Product datasheet for **TA368592S**


KANK3 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human KANK3
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	KN motif and ankyrin repeat domains 3
Database Link:	Entrez Gene 256949 Human Q6NY19
Background:	May be involved in the control of cytoskeleton formation by regulating actin polymerization.
Synonyms:	ANKRD47; FLJ46061



[View online »](#)

Product images:


Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA368592] (KANK3 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA368592] (KANK3 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: ×200)