

Product datasheet for **TA369059S**

Cyclophilin B (PPIB) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human fetal liver tissue and 293T cells, human liver cancer tissue and hela cells IHC: 100-300 Positive control: Human cervical cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human PPIB
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	24 kDa
Gene Name:	peptidylprolyl isomerase B
Database Link:	Entrez Gene 5479 Human P23284
Background:	The protein encoded by this gene is a cyclosporine-binding protein and is mainly located within the endoplasmic reticulum. It is associated with the secretory pathway and released in biological fluids. This protein can bind to cells derived from T- and B-lymphocytes, and may regulate cyclosporine A-mediated immunosuppression. Variants have been identified in this protein that give rise to recessive forms of osteogenesis imperfecta.
Synonyms:	CYP-S1; CYPB; MGC2224; MGC14109; OI9; PPIase; rotamase; S-cyclophilin; SCYLP




[View online »](#)


Product images:




Gel: 10%SDS-PAGE
Lysate: 40 µg
Lane 1-4: Human fetal liver tissue
293T cells
human liver cancer tissue
hela cells
Primary antibody: [TA369059] (PPIB Antibody) at
dilution 1/450
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 20 seconds




Immunohistochemistry of paraffin-embedded
Human cervical cancer tissue using [TA369059]
(PPIB Antibody) at dilution 1/40 (Original
magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA369059] (PPIB Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369059] (PPIB Antibody) at dilution 1/40 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369059] (PPIB Antibody) at dilution 1/40, treated with fusion protein. (Original magnification: $\times 200$)