

Product datasheet for **TA369198S**

LZTFL1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human testis tissue, 293T and A 172 cells IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human LZTFL1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	35 kDa
Gene Name:	leucine zipper transcription factor like 1
Database Link:	Entrez Gene 54585 Human Q9NQ48




[View online »](#)

Background:

This gene encodes a ubiquitously expressed protein that localizes to the cytoplasm. This protein interacts with Bardet-Biedl Syndrome (BBS) proteins and, through its interaction with BBS protein complexes, regulates protein trafficking to the ciliary membrane. Nonsense mutations in this gene cause a form of Bardet-Biedl Syndrome; a ciliopathy characterized in part by polydactyly, obesity, cognitive impairment, hypogonadism, and kidney failure. This gene may also function as a tumor suppressor; possibly by interacting with E-cadherin and the actin cytoskeleton and thereby regulating the transition of epithelial cells to mesenchymal cells. Alternative splicing of this gene results in multiple transcript variants.

Synonyms:

FLJ36386

Product images:

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane 1-3: Human testis tissue


293T cells

A172 cells


Primary antibody: [TA369198] (LZTFL1 Antibody)
at dilution 1/300Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 20 seconds


Immunohistochemistry of paraffin-embedded
Human liver cancer tissue using [TA369198]
(LZTFL1 Antibody) at dilution 1/30 (Original
magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369198] (LZTFL1 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA369198] (LZTFL1 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA369198] (LZTFL1 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)