

## Product datasheet for **TA369363S**


### SNAP47 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 200-400 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human SNAP47
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	synaptosome associated protein 47kDa
Database Link:	<a href="#">Entrez Gene 116841 Human</a> <a href="#">Q5SQN1</a>
Background:	Plays a role in intracellular membrane fusion.
Synonyms:	C1orf142; DKFZp686M10160; FLJ12517; HEL170; SNAP-47; SVAP1



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA369363] (SNAP47 Antibody) at dilution 1/200 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA369363] (SNAP47 Antibody) at dilution 1/200, treated with fusion protein. (Original magnification:  $\times 200$ )