

## Product datasheet for **TA369399S**


### THAP10 Rabbit Polyclonal Antibody

#### Product data:


|                       |                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                                                                                                                                                                                                                                                                                       |
| Applications:         | IHC                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilution: | IHC: 25-100<br>Positive control: Human liver cancer<br>Predicted cell location: Cytoplasm                                                                                                                                                                                                                                                                |
| Reactivity:           | Human                                                                                                                                                                                                                                                                                                                                                    |
| Host:                 | Rabbit                                                                                                                                                                                                                                                                                                                                                   |
| Isotype:              | IgG                                                                                                                                                                                                                                                                                                                                                      |
| Clonality:            | Polyclonal                                                                                                                                                                                                                                                                                                                                               |
| Immunogen:            | Full length fusion protein                                                                                                                                                                                                                                                                                                                               |
| Formulation:          | pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol                                                                                                                                                                                                                                                                                                         |
| Purification:         | Antigen affinity purification                                                                                                                                                                                                                                                                                                                            |
| Conjugation:          | Unconjugated                                                                                                                                                                                                                                                                                                                                             |
| Storage:              | Store at -20°C.                                                                                                                                                                                                                                                                                                                                          |
| Stability:            | 1 year                                                                                                                                                                                                                                                                                                                                                   |
| Gene Name:            | THAP domain containing 10                                                                                                                                                                                                                                                                                                                                |
| Database Link:        | <a href="#">Entrez Gene 56906 Human Q9P2Z0</a>                                                                                                                                                                                                                                                                                                           |
| Background:           | This gene encodes a member of a family of proteins sharing an N-terminal Thanatos-associated domain. The Thanatos-associated domain contains a zinc finger signature similar to DNA-binding domains. This gene is part of a bidirectional gene pair on the long arm of chromosome 15 that is regulated by estrogen and may play a role in breast cancer. |
| Synonyms:             | THAP10                                                                                                                                                                                                                                                                                                                                                   |



[View online »](#)

**Product images:**


Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369399] (THAP10 Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369399] (THAP10 Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)