

Product datasheet for **TA369479S**


XAB2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human prostate cancer Predicted cell location: Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human XAB2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	XPA binding protein 2
Database Link:	Entrez Gene 56949 Human Q9HCS7
Background:	Involved in transcription-coupled repair (TCR), transcription and pre-mRNA splicing.
Synonyms:	DKFZp762C1015; HCNP; HCRN; KIAA1177; NTC90; SYF1



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA369479] (XAB2 Antibody) at dilution 1/25 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA369479] (XAB2 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: $\times 200$)