

Product datasheet for **TA369630S**

LYPD4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 10-50 Positive control: Human prostate cancer Predicted cell location: Cell membrane
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	LY6/PLAUR domain containing 4
Database Link:	Entrez Gene 147719 Human Q6UWN0
Background:	The specific function of LYPD4 is not yet known.
Synonyms:	MGC42718; SMR




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA369630] (LYPD4 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA369630] (LYPD4 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)