

## Product datasheet for **TA369654S**

### **ZNF313 (RNF114) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-200 Positive control: Human thyroid cancer Predicted cell location: Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human RNF114
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	ring finger protein 114
Database Link:	<a href="#">Entrez Gene 55905 Human Q9Y508</a>




[View online »](#)

**Background:**


E3 ubiquitin-protein ligase that promotes the ubiquitination of various substrates (PubMed:23645206, PubMed:25165885). In turn, participates in the regulation of many biological processes including cell cycle, apoptosis, osteoclastogenesis as well as innate or adaptive immunity (PubMed:25165885, PubMed:28708287). Acts as negative regulator of NF-kappa-B-dependent transcription by promoting the ubiquitination and stabilization of the NF-kappa-B inhibitor TNFAIP3 (PubMed:25165885). May promote the ubiquitination of TRAF6 as well (PubMed:28708287). Acts also as a negative regulator of T-cell activation (PubMed:25165885). Inhibits cellular dsRNA responses and interferon production by targeting MAVS component for proteasomal degradation (PubMed:25165885). Ubiquitinates the CDK inhibitor CDKN1A leading to its degradation and probably also CDKN1B and CDKN1C (PubMed:23645206). This activity stimulates cell cycle G1-to-S phase transition and suppresses cellular senescence. May play a role in spermatogenesis.

**Synonyms:**

ZNF228; ZNF313

**Product images:**

Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA369654] (RNF114 Antibody) at dilution 1/90 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA369654] (RNF114 Antibody) at dilution 1/90, treated with fusion protein. (Original magnification: ×200)